

Laboratory Name _____

Address _____

City _____ State _____ ZIP _____

Phone _____

E-mail _____

Patient Name _____ Case # _____





NOTE: CMC offers free next-day delivery on all abutment cases. Lab screws included.

ABUTMENT MATERIAL CHOICE (please check box for the material of your choice for each tooth #)

Titanium	Gold-Hue Anodized	Zirconia	Chrome Cobalt	Full-Contour Screw-Retained	Tooth #	Implant Brand	Platform Φ
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	_____
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	_____
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	_____
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	_____
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	_____
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	_____

EMERGENCE WIDTH OPTION (select one)

- Full anatomical dimensions
 Contour soft tissue (default)
 Support tissue
 No tissue displacement



MARGINS

Buccal / facial : _____
 Distal : _____
 Mesial : _____
 Lingual : _____

DEFAULT

1.0 mm subgingival
 0.75 mm subgingival
 0.75 mm subgingival
 0.5 mm subgingival

CLINICIAN SPECIFIED (if different from default)

Buccal / facial : _____
 Distal : _____
 Mesial : _____
 Lingual : _____

SPECIAL INSTRUCTIONS

Screen shot for design approval needed? YES NO

DUE DATE _____ / _____ / _____

PROMO CODE _____

SIGNATURE _____ **DATE** _____

For zirconia abutments: All internal connections consist of a two piece hybrid system with titanium base. All external connections consist of a one piece zirconia system.